

Eletrodos PadPro ConMed

Eletrodo de Desfibrilação / Marca passo / Cardioversão / Monitoração

- Melhor posicionamento ao corpo do paciente
- Excelente aderência em pacientes com pêlos
- Distribuição uniforme de energia em toda a superfície do eletrodo
- Eliminação de possíveis "tenting" causador de faíscas e queimaduras
- Tamanho do eletrodo adulto: 15,2 cm x 10,8 cm
- Tamanho do eletrodo pediátrico: 10,8 cm x 7,3 cm
- Validade (embalagem lacrada) 24 meses a contar da data de fabricação
- Compatível com desfibriladores monofásicos e bifásicos
- Livre de látex
- Fabricado nos USA por CONMED



Adulto



Pediátrico

Conectores:



Philips



Philips



Zoll (M series)








Universal



Medtronic
Physio Control

Tabela de compatibilidade

 <p>Anderson Connector ¹</p>		<p>1- Electrodes with the Anderson connector can be utilized in environments with multiple brands of defibrillators (via PadPro® adapters or PadPro® cables). Cmos-Drake / Teb CV-10MP</p>
 <p>Zoll® ²</p>		<p>2 - Zoll M-Series models are comparable to Zoll® multifunction electrodes.</p>
 <p>Medtronic® ³</p>		<p>3 - Medtronic models are comparable to Medtronic Physio Control® QUIK-COMBO™ electrodes.</p>
 <p>Philips Barrel-Style ⁴</p>		<p>4 - Philips Barrel-style models are comparable to Philips Codemaster (Barrel-Style) electrodes.</p>
 <p>Philips Plug-Style ⁵</p>		<p>5 - Philips Plug-Style models are comparable to Philips Heartstart (Plug-Style) electrodes.</p>

* Product and company names listed are trademarks or trade names of their respective companies.